



Bandage Ontario 2088

Indication


- After thorax surgery
- After abdominal surgery
- For thorax trauma such as thoracic contusions or broken ribs
- For sciatic problems
- After disc herniation
- After weight loss in the case of obesity

Application

- Regular compression distribution
- Helps the tummy regain its shape after birth
- Breathing is unobstructed, even with a tight fit
- Flexibly adjustable during the healing process

Breathable, softly padded wound dressing covering

Height: 25 cm



Steplessly adjustable using
hook & loop fastener

Simple and quick application/
changing of the wound dressing

Width and length adjustment: 10 cm

SIZE CIRCUMFERENCE WAIST

0	50 – 65 cm
1	65 – 80 cm
2	80 – 95 cm
3	95 – 105 cm
4	105 – 115 cm
5	115 – 125 cm
6	125 – 140 cm
7	140 – 150 cm

AVAILABLE COLOURS

006 white



FABRIC

55% PA, 16% PES, 13% CO, 8% EL, 8% CV

95° C washable